


MP-30T Syringe Driver

The MP-30T Syringe Driver is **small, smart** and most importantly **safe**

The MP-30T Syringe Driver is compact, lightweight and space-saving. Featuring touch-screen technology, with visual and audible alarms, the MP-30T is user friendly and benefits from adjustable occlusion levels and real-time monitoring of infusion states.



- 1 Adjustable occlusion levels
- 2 Visual and audible alarm information
- 3 Convenient touch screen with clear interface
- 4 Stackable between any units in the MEDCAPTAIN MP range (11-5985 & 11-5990)
- 5 User defined screen brightness and background colours
- 6 2 year warranty



Pumps in the MP range can be connected and secured easily



A slide-on handle is supplied with every pump for portability



A clamp is supplied with each pump to connect to an Infusion Stand or Kennel

Specifications:

Dimensions:	244 (W) x 74 (H) x 164 (D) mm
Weight:	<1.6kg (including battery)
Screen:	2.4" LCD touch-screen, 7 colours selectable
Drug library:	Up to 2000 drugs can be built in
Power supply switching:	Automatically switches to battery when power supply is cut off
Compatibility:	All syringes: 5 ml, 10 ml, 20 ml, 30 ml, 50 ml/60 ml standard
Battery:	Lithium 7.2v, 2600mAh
Battery duration:	> 6 hours (based on 5 ml/h with a new battery)

Please Note: WIFI option available, please call our sales team for more information



Basic Parameter:

Infusion mode	Rate, Time, Weight, Sequence, Relay mode, Loading dose, trapezia mode, micro mode, Tiva	
Infusion rate range	0.1-100.0 ml/h (5ml syringe) 0.1-600.0 ml/h (20ml syringe) 0.1-2000.0 ml/h (50ml syringe)	0.1-300.0 ml/h (10ml syringe) 0.1-900.0 ml/h (30ml syringe)
Increment	0.1 – 99.99 ml/h, (min. increment 0.01 ml/h) 100– 999.9 ml/h, min. increment 0.1 ml/h 1000-1200 ml/h, min. increment 1 ml/h	
VTBI range	0.1 – 99.99 ml (minimum increment 0.01) 100 – 999.9 ml (minimum increment 0.1) 1000 – 9999 ml (minimum increment 1)	
Pre-set time	00.00 – 99.59 (hh:mm) adjustable , step by 1 min.	
Accumulated volume	0.0 – 99999.99 ml	
Accuracy	Mechanical accuracy: ±1%, Accuracy including syringe: ±2%	
KVO rate	.01 – 5.0 ml/h adjustable, step by 0.01 ml/h, default rate: 1 ml/h	
Bolus rate	0.1-100.0 ml/h (5 ml syringe) 0.1-300.0 ml/h (10 ml syringe) 0.1-600.0 ml/h (20 ml syringe) 0.1-900.0 ml/h (30 ml syringe) 0.1-2000.0 ml/h (50 ml/60 ml syringe) 3 bolus ways selectable: Automatic/ Manual/ Rapid quantitative Bolus	
Purge range	100.0 ml/h (5 ml syringe) 300.0 ml/h (10 ml syringe) 600.0 ml/h (20 ml syringe) 900.0 ml/h (30 ml syringe) 2000.0 ml/h (50 ml/60 ml syringe)	
Occlusion level	11 level selectable, 225 to 975mmHg	

01347 878697 | sales@jakmarketing.co.uk
www.jakmarketing.co.uk

J.A.K Marketing Limited Perfect House, 11-21 Dale Road, Sheriff Hutton, York, YO60 6RZ

Registered Company No. 4316591

